



Ultra-Fast GC for Vocus PTR- TOF

Separation of isomeric species and enhanced chemical specificity for complex mixtures



Features

- Compact design with plug-and-play coupling to Vocus PTR-TOF
- Detect trace species with thermal desorption preconcentrator trap
- Results in as fast as 2 min
- Fast switching between continuous measurements and GC separations
- Integrated software control with customizable methods and workflows
- Large variety of column and traps targeting specific application

Applications

- Air monitoring
- Food, flavor, and fragrances
- Beverage and alcohol
- Material emissions
- Forensics
- Petrochemical analysis

TOFWERK
www.tofwerk.com
ptr@tofwerk.com

